



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bookwalter Urology Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Bookwalter Urology Retractor System

SKU RM-BURS-1

\$2,399 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Field Post	12" (30.5 cm)
1	Post Coupler	Standard
1	Horizontal Flex Bar	1" Diameter (2.5 cm)
8	Tilt Ratchet Mechanisms	Multi-Angle
1	Segmented Ring	10½" × 19" (26.7 × 48.3 cm)
1	Segmented Ring	8½" × 14" (21.6 × 35.6 cm)
1	Malleable Blade	2" × 6" (5.1 × 15.2 cm)
1	Malleable Blade	1½" × 6" (3.8 × 15.2 cm)
2	Kelly Blades	1" × 2½" (2.5 × 6.4 cm)
1	Kelly Blade	1" × 1½" (2.5 × 3.8 cm)
2	Kelly Blades	2" × 2½" (5.1 × 6.4 cm)
3	Vaginal Lateral Blades	1" × 3" (2.5 × 7.6 cm)
1	Renal Wylie Blade	1" × 7" (2.5 × 17.8 cm)
1	Sterilization Box	Protective Case

PRODUCT DESCRIPTION

The Bookwalter Urology Retractor System is a table-mounted self-retaining retractor set designed to provide stable, hands-free exposure during abdominal and flank approaches to the kidneys, ureters, bladder and retroperitoneal urological structures. The Bookwalter Urology Retractor System mounts to the operating table with a 12" (30.5 cm) field post, a standard post coupler and a 1" (2.5 cm) diameter horizontal flex bar.

Two segmented rings are supplied — 10½" × 19" (26.7 × 48.3 cm) and 8½" × 14" (21.6 × 35.6 cm) — so the frame can be matched to the incision and the patient. Eight multi-angle tilt ratchet mechanisms secure the blades to the ring and allow retraction depth and blade angle to be set independently.

Retraction hardware covers deep and superficial exposures: 2" × 6" and 1½" × 6" malleable blades, Kelly blades in 1" × 2½", 1" × 1½" and 2" × 2½" sizes, three 1" × 3" vaginal lateral blades and a 1" × 7" renal Wylie blade. All components are manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and supplied in a sterilization box.

The Bookwalter Urology Retractor System is supplied non-sterile, is reusable and fully autoclavable, and must be cleaned and sterilized according to your facility's validated procedures before use. Retractor Maker manufactures and supplies the set directly from Sialkot, Pakistan.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)