



RETRACTOR MAKER

PRODUCT BROCHURE OF

Lateral Ebonized Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Lateral Ebonized Vaginal Retractor

SKU RM-LEVR-63-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Lateral Ebonized Vaginal Retractor	Dimensions / Details: 6.3 cm size parameter, laser-safe ebonized blade · Quantity Included: 1	\$79
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$79

PRODUCT DESCRIPTION

The Lateral Ebonized Vaginal Retractor is a specialized gynecological surgical instrument engineered specifically to provide safe lateral vaginal wall retraction during laser surgery and electrosurgical procedures. Coated with a non-reflective black ebonized finish, this retractor absorbs stray laser radiation and eliminates laser beam reflection, protecting surrounding healthy vaginal mucosa from accidental thermal injury.

Built for procedural safety and durability, this Lateral Ebonized Vaginal Retractor configuration features a standard surgical size parameter of 6.3 cm (\$79). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its surface profile finish combines highly polished handling interfaces with a laser-safe ebonized working surface layer that withstands high-temperature ablation environments. The curved blade contours displace the lateral vaginal walls smoothly, providing an unhindered visual corridor to the cervix and fornices.

Manufactured to strict clinical tolerances, Retractor Maker produces the Lateral Ebonized Vaginal Retractor to withstand repetitive hospital autoclave cycles, delivering exceptional procedural hygiene, high corrosion resistance, and complete clinical durability for laser gynecological suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)