



RETRACTOR MAKER

PRODUCT BROCHURE OF

Senn Green Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Senn Green Retractor

SKU RM-SGR-14-1

\$27 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
6 mm × 10 mm Blade — 14 cm	RM-SGR-14-1-1	\$27
6 mm × 20 mm Blade — 14 cm	RM-SGR-14-1-2	\$29

PRODUCT DESCRIPTION

The Senn Green Retractor is a specialized handheld surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and strap muscle retraction during thyroid, endocrine, and superficial plastic surgery procedures. Combining an elongated right-angled blade with a downward-curved retaining lip, this instrument allows surgical assistants to maintain constant counter-traction on delicate tissue planes without causing local bruising or lacerations.

Built to withstand frequent autoclave sterilization cycles while remaining highly resistant to corrosion, this Senn Green Retractor maintains an overall tool length of 14 cm. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument is highly adaptable to varying incision depths through two distinct blade footprint options selectable via the menu parameters: a 6 mm × 10 mm blade layout (\$27) for superficial wound margins, and an extended 6 mm × 20 mm blade layout (\$29) for deeper cervical and muscular exposure. Its flat ergonomic handle ensures fatigue-free manual grip stability.

Classified under certified premium operating room grade hooks and handheld retractors, Retractor Maker produces the Senn Green Retractor with an anti-glare satin matte finish, delivering exceptional structural rigidity, corrosion resistance, and complete autoclave durability for hospital operating theaters and outpatient surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)