



RETRACTOR MAKER

PRODUCT BROCHURE OF

Pediatric Ribbon Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Pediatric Ribbon Retractor

SKU RM-PRR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
17.5 cm (6.9")	RM-PRR-1-1	\$44
20.0 cm (7.9")	RM-PRR-1-2	\$49
22.0 cm (8.7")	RM-PRR-1-3	\$54

PRODUCT DESCRIPTION

The Pediatric Ribbon Retractor is an indispensable malleable surgical retractor engineered specifically to provide safe, non-traumatic tissue displacement and visceral shielding in neonatal and pediatric operative fields. Designed with a highly ductile flat ribbon blade, this instrument allows pediatric surgeons to shape and bend the spatula intraoperatively to match unique infant anatomical curves, avoiding localized pressure trauma on delicate viscera.

Built with an adaptive selection matrix combining customized tool lengths and blade widths, this Pediatric Ribbon Retractor series is available in three overall lengths: 17.5 cm (\$44), 20.0 cm (\$49), and 22.0 cm (\$54). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), surgeons can select across three working widths: 8 mm for narrow corridors, 10 mm for standard infant procedures, and 12 mm for broader abdominal wall holding. The smoothly polished margins prevent serosal tears, while the malleable body retains its formed shape under tissue counter-traction without springback.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Pediatric Ribbon Retractor to withstand repetitive intraoperative bending and hospital sterilization cycles, providing pediatric surgery centers and congenital heart suites worldwide with exceptional ductility, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)