



RETRACTOR MAKER

PRODUCT BROCHURE OF

Fischl Shoe Horn Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Fischl Shoe Horn Retractor

SKU RM-FSHR-16-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Fischl Shoe Horn Retractor	16 cm overall length, curved shoe-horn spatula blade	\$24
Protective Sleeve / Box	Individual transit packaging	\$24

PRODUCT DESCRIPTION

The Fischl Shoe Horn Retractor is a specialized handheld surgical instrument engineered specifically to provide smooth, atraumatic soft-tissue displacement and tissue plane separation across plastic, general, and reconstructive procedures. Designed with an ultra-smooth, contoured shoe-horn spatula blade, this instrument allows surgeons to guide implants, elevate delicate skin flaps, and retract wound margins without causing friction trauma or edge contusions.

Built for versatile clinical utility, this Fischl Shoe Horn Retractor configuration features an overall length of 16 cm (\$24). Its working blade incorporates a gently scooped, curved profile that slips effortlessly into narrow surgical tunnels, providing broad-surface counter-traction that prevents focal pressure necrosis. The ergonomic flat handle delivers exceptional manual tactile feedback, allowing operating assistants to maintain stable tissue deflection throughout extended dissections.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Fischl Shoe Horn Retractor to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for hospital operating theaters and cosmetic surgery clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)