



RETRACTOR MAKER

PRODUCT BROCHURE OF

Taylor Spinal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Taylor Spinal Retractor

SKU RM-TSR-1

\$74 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
7 ¼" (18.42 cm) — 3" (7.62 cm) depth — Cervical and thoracic posterior laminectomy, adolescent spine	RM-TSR-1-1	\$74
7 ½" (19.05 cm) — 3" or 4" depth — Standard adult lumbar spine laminectomy and discectomy	RM-TSR-1-2	\$79
8.0" (20.32 cm) — 4" (10.16 cm) depth — Deep lumbar fusion, bariatric spine, and sacral exposures	RM-TSR-1-3	\$89

PRODUCT DESCRIPTION

The Taylor Spinal Retractor is a classic spinal surgical retractor engineered specifically to provide powerful, non-slip paraspinal muscle displacement and spinous process anchoring during posterior spinal procedures. Designed with an elongated right-angled blade terminating in a sharp spike tip, this instrument hooks directly over the spinous process or anchors into intervertebral bone margins, holding thick paraspinal musculature laterally away from the lamina.

Built to withstand frequent autoclave sterilization cycles while remaining highly resistant to corrosion, this Taylor Spinal Retractor series is offered in three overall lengths: 7 ¼ inches (18.42 cm, \$74), 7 ½ inches (19.05 cm, \$79), and 8 inches (20.32 cm, \$89). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), surgeons can select between two right-angled blade sizes: a 3-inch (7.62 cm) reach depth for standard lumbar corridors or a 4-inch (10.16 cm) reach depth for deep spinal canals and muscular patients. Its handle assembly layout features a single-ended design integrated into a solid L-shaped handle configuration built for optimal palm placement, high traction control, and towel weight attachment.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Taylor Spinal Retractor to withstand immense manual pulling loads, delivering exceptional structural rigidity, corrosion resistance, and complete autoclave durability for spine centers, neurosurgical suites, and orthopedic operating rooms worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)