



RETRACTOR MAKER

PRODUCT BROCHURE OF

Low Profile Foot and Ankle Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Low Profile Foot and Ankle Retractor

SKU RM-LPFAR-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
7.5 inches	RM-LPFAR-1-1	\$84
8.75 inches	RM-LPFAR-1-2	\$89

PRODUCT DESCRIPTION

The Low Profile Foot and Ankle Retractor is an indispensable orthopedic surgical instrument engineered specifically to provide dependable bone leverage and soft-tissue displacement in tight podiatric and ankle joint spaces. Designed with a low-profile shank, this instrument allows surgeons to work comfortably in confined operative fields without bulky retractor bodies obstructing surgical drills, screws, or fluoroscopic imaging.

Built to accommodate varying tissue depths, the Low Profile Foot and Ankle Retractor series is offered in two overall lengths: 7.5 inches (19.05 cm, \$84) for forefoot and midfoot procedures, and 8.75 inches (22.22 cm, \$89) for deeper hindfoot and ankle joint exposures. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers structural blade profiles in both straight and bent (curved) geometry layouts to adapt to complex osseous borders. Its contoured beak tip seats securely behind cortical bone, while the ergonomic handle delivers steady manual counter-traction.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Low Profile Foot and Ankle Retractor to withstand heavy orthopedic torque and repeated hospital autoclave cycles, providing podiatric, trauma, and orthopedic foot and ankle specialists worldwide with unyielding clinical durability and surgical precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)