



RETRACTOR MAKER

PRODUCT BROCHURE OF

Meyerding Skin Hook/Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Meyerding Skin Hook/Retractor

SKU RM-MSHR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
16 cm (approx. 6.3")	RM-MSHR-1-1	\$44
22.86 cm (9.0")	RM-MSHR-1-2	\$49

PRODUCT DESCRIPTION

The Meyerding Skin Hook/Retractor is a versatile double-ended surgical instrument engineered specifically to combine fine cutaneous traction with deeper soft-tissue displacement in a single chassis. Designed for high procedural agility, this dual-ended tool allows surgical teams to elevate skin edges cleanly with a fine hook on one end, then reverse the instrument to hold deeper fascial margins with a solid retractor blade without exchanging tools.

Built for complete procedural adaptability, this Meyerding Skin Hook/Retractor configuration is offered in two overall instrument lengths: a compact 16 cm model (\$44) for fine cosmetic procedures and an extended 22.86 cm model (\$49) for deeper operative fields. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the integrated working ends feature a Retractor End equipped with a solid 5.0 mm wide × 16.0 mm long blade profile, paired with a Skin Hook End featuring 1 sharp prong to ensure precise, slippage-free engagement with sensitive tissue margins. Its balanced central handle provides fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Meyerding Skin Hook/Retractor to withstand repetitive hospital sterilization cycles, providing plastic, reconstructive, and general surgical suites worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)