



RETRACTOR MAKER

PRODUCT BROCHURE OF

Beckmann Adson Self-Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Beckmann Adson Self-Retaining Retractor

SKU RM-BASRR-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
25 mm Shaft Profile	RM-BASRR-1-1	\$49
44 mm Shaft Profile	RM-BASRR-1-2	\$59

PRODUCT DESCRIPTION

The Beckmann Adson Self-Retaining Retractor is a heavy-duty surgical retractor engineered specifically to provide broad-surface, hands-free muscular retraction during deep spinal and suboccipital procedures. Designed with an elongated 31 cm frame, this instrument allows surgeons to maintain deep working corridors while keeping the ratchet handles safely outside the direct line of surgical vision.

Built for maximum tissue engagement under heavy mechanical tension, this Beckmann Adson Self-Retaining Retractor model features an expansive 5×4 prong array that distributes traction forces evenly across dense muscle planes to prevent localized tearing. Surgeons can select between sharp prongs for firm fascial grip or blunt prongs to protect neurovascular tissue, paired with either a 25 mm (\$49) or 44 mm (\$59) shaft diameter profile. Its finger-ring ratchet mechanism locks securely at micro-incremental steps.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Beckmann Adson Self-Retaining Retractor to withstand heavy counter-traction, providing spinal, orthopedic, and neurosurgical departments with reliable durability and full autoclave compatibility.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)