



RETRACTOR MAKER

PRODUCT BROCHURE OF

Saldanha-Solz Abdominal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Saldanha-Solz Abdominal Retractor

SKU RM-SSAR-235-1

\$199 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Saldanha-Solz Abdominal Retractor	23.5 cm overall length, tissue-optimized curved blade	\$199
Protective Sleeve / Box	Individual transit packaging	\$199

PRODUCT DESCRIPTION

The Saldanha-Solz Abdominal Retractor is a specialized aesthetic and reconstructive plastic surgery retractor engineered specifically to provide expansive soft-tissue elevation and superior visual clearance during lipo-abdominoplasty and high-definition body sculpting. Designed to support thick abdominal cutaneous flaps, this instrument allows plastic surgeons to maintain wide operative pockets without crushing dermal vascularity or tearing fragile adipose lobules.

Built to withstand frequent autoclave sterilization cycles while providing a lifetime warranty against rust and defects, this Saldanha-Solz Abdominal Retractor configuration features an overall tool length of 23.5 cm (\$199). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry layout features a single-ended, tissue-optimized retractor blade contoured explicitly to protect abdominal structures during fat harvesting and muscle plication maneuvers. Its heavy-duty ergonomic handle delivers fatigue-free manual leverage for operating assistants throughout prolonged body contouring cases.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Saldanha-Solz Abdominal Retractor to withstand repetitive hospital sterilization cycles, providing cosmetic, plastic, and reconstructive surgery centers worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)