



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bookwalter General Surgery Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Bookwalter General Surgery Retractor

SKU RM-BGSR-1

\$1,999 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Field Post	12" (30.5 cm)
1	Horizontal Bar	Standard
1	Post Coupler	Standard
1	Horizontal Flex Bar	1" Diameter (2.5 cm)
2	Standard Ratchets	Non-Tilt
6	Tilt Ratchet Mechanisms	Multi-Angle
1	Oval Ring, Medium	10½" × 14½" (26.7 × 36.8 cm)
1	Segmented Ring, Medium	10½" × 19" (26.7 × 48.3 cm)
1	Balfour Blade	2⅝" × 3" (6.7 × 7.6 cm)
2	Kelly Blades	2" × 3" (5.1 × 7.6 cm)
1	Kelly Blade	2" × 4" (5.1 × 10.2 cm)
1	Kelly Blade	2" × 5" (5.1 × 12.7 cm)
1	Kelly Blade	2" × 6" (5.1 × 15.2 cm)
1	Malleable Blade	3" × 6" (7.6 × 15.2 cm)

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
2	Malleable Blades	2" × 6" (5.1 × 15.2 cm)
1	Malleable Blade	1½" × 6" (3.8 × 15.2 cm)
1	Harrington Blade	1" × 5⅝" × 2½" (2.5 × 14.3 × 6.4 cm)
1	Gelpi Blade	Standard
1	Sterilization Box	Protective Case

PRODUCT DESCRIPTION

The Bookwalter General Surgery Retractor is a table-mounted self-retaining retractor designed to provide stable, hands-free exposure during open abdominal and general surgical procedures. Its modular framework combines a field post, horizontal bar, horizontal flex bar, post coupler, retractor rings and adjustable ratchet mechanisms to support selected blades around the operative field.

The General Surgery configuration includes a medium oval ring, a medium segmented ring, two standard non-tilt ratchets and six multi-angle tilt ratchet mechanisms. The set also includes Balfour, Kelly, malleable, Harrington and Gelpi blades in different configurations.

The complete system is supplied in a dedicated sterilization box for organized storage and processing. The Bookwalter General Surgery Retractor is supplied non-sterile and is designed to be reusable and autoclavable. Cleaning and sterilization should be performed according to the healthcare facility's validated procedures before use.

Retractor Maker supplies the Bookwalter General Surgery Retractor directly to hospitals, distributors and other B2B customers.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)