



RETRACTOR MAKER

PRODUCT BROCHURE OF

Children's Hospital Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Children's Hospital Retractor

SKU RM-CHR-PED-1

\$39 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Children's Hospital Retractor	14 cm overall length, double-ended pediatric blades	\$39
Protective Sleeve / Box	Individual transit packaging	\$39

PRODUCT DESCRIPTION

The Children's Hospital Retractor is a specialized pediatric surgical retractor engineered specifically to provide gentle, atraumatic tissue retraction in delicate neonatal, pediatric, and superficial adult procedures. Designed with a functional double-ended architecture, this instrument gives pediatric surgeons two distinct blade options in a single tool, optimizing operative efficiency and reducing surgical tray clutter.

Built for precision in localized operative fields, this Children's Hospital Retractor model features an overall length of 14 cm (\$39). Its working ends incorporate smooth, polished blades with rounded edges specifically contoured to prevent focal pressure necrosis and accidental tissue tearing in fragile pediatric anatomies. The lightweight, ergonomically balanced central handle allows surgical assistants to maintain steady counter-traction with minimal fatigue.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Children's Hospital Retractor to withstand repetitive hospital sterilization cycles, providing pediatric surgery centers, children's hospitals, and general operating theaters with high corrosion resistance and clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)