



RETRACTOR MAKER

PRODUCT BROCHURE OF

Maxwell Flap Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Maxwell Flap Retractor

SKU RM-MFR-4IN-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
22 mm wide layout	RM-MFR-4IN-1-1	\$29
38 mm wide layout	RM-MFR-4IN-1-2	\$34

PRODUCT DESCRIPTION

The Maxwell Flap Retractor is a specialized handheld plastic and reconstructive surgery retractor engineered specifically to provide firm, non-slip skin flap elevation during rhytidectomy, blepharoplasty, and soft-tissue procedures. Designed with sharp rake-style prongs, this instrument allows surgeons to apply uniform manual counter-traction across cutaneous and subcutaneous planes without causing focal dermal tearing.

Built for delicate manual control in tight operative spaces, this Maxwell Flap Retractor maintains a uniform body length of 4 inches (10.2 cm) that provides operating assistants with exceptional fingertip sensitivity. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working blade incorporates a 4 sharp rake-style prongs layout available across two selectable working end widths: a 22 mm wide layout (\$29) for localized facial flaps and an expanded 38 mm wide layout (\$34) for broad cheek and neck undermining. Its ergonomic flat handle ensures fatigue-free grip stability.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Maxwell Flap Retractor to withstand repetitive hospital autoclave cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for cosmetic and facial plastic surgery clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)