



RETRACTOR MAKER

PRODUCT BROCHURE OF

Breast Engler Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Breast Engler Retractor

SKU RM-BER-OPT-1

\$139 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Breast Engler Retractor	65 × 12 mm blade with integrated fiber-optic bundle
1	Universal Light Cable Fitting	Connects to standard OR light cables
1	Protective Sleeve / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Breast Engler Retractor is an indispensable illuminated surgical retractor engineered specifically to provide superior pocket visualization and gentle parenchymal elevation during aesthetic and reconstructive breast surgery. Designed with an integrated fiber-optic bundle, this instrument projects cool, shadow-free illumination directly into deep subglandular and subpectoral planes, eliminating dark recesses without requiring clumsy external headlamps.

Built to match delicate breast anatomy, this Breast Engler Retractor configuration features a precision-crafted 65 mm long by 12 mm wide blade profile (\$139). Its narrow, elongated architecture slips effortlessly through minimal inframammary or periareolar incisions while providing unyielding mechanical leverage to lift pectoral muscle flaps. The universal fiber-optic connector integrates seamlessly with standard operating room light source cables, ensuring versatile clinical connectivity.

Crafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-transmission optical glass fibers with an anti-glare satin finish, Retractor Maker supplies the Breast Engler Retractor non-sterile and fully autoclavable inside protective packaging, providing cosmetic, plastic, and oncologic breast surgeons with high durability and surgical clarity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)