



RETRACTOR MAKER

PRODUCT BROCHURE OF

Combination Wound and Cheek Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Combination Wound and Cheek Retractor

SKU RM-CWCR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Combination Wound and Cheek Retractor	Double-ended oral and wound retractor — 1	\$69
Protective Sleeve / Box	Individual transit packaging — 1	\$69

PRODUCT DESCRIPTION

The Combination Wound and Cheek Retractor is an adaptable, dual-purpose surgical instrument engineered to deliver superior soft-tissue displacement during oral, maxillofacial, and plastic reconstructive procedures. Designed with contrasting functional ends, this retractor allows surgical teams to transition effortlessly between broad cheek elevation and localized wound margin retraction without changing instruments.

Built for optimal intraoral and facial clearance, this Combination Wound and Cheek Retractor configuration (\$69) integrates a broad, smoothly contoured cheek blade on one side to protect the buccal mucosa from rotating burrs, and a curved hook blade on the opposing end for precise skin flap and wound retraction. Its balanced central handle provides surgical assistants with steady manual control and comfortable grip ergonomics, reducing hand fatigue during lengthy procedures.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the Combination Wound and Cheek Retractor non-sterile and fully autoclavable inside protective packaging, delivering high corrosion resistance and dependable durability for operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)