



RETRACTOR MAKER

PRODUCT BROCHURE OF

Dockhorn Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Dockhorn Retractor

SKU RM-DOCK-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
12.0 mm × 40.0 mm blade	RM-DOCK-1-1	\$64
14.0 mm × 55.0 mm blade	RM-DOCK-1-2	\$69

PRODUCT DESCRIPTION

The Dockhorn Retractor is a specialized handheld surgical retractor engineered specifically to provide accurate, stable soft-tissue displacement and wound margin isolation across oral, maxillofacial, and general surgical procedures. Designed with an ergonomic shank and a distinctive concave upper blade surface, this instrument holds wound margins firmly while directing fluid and reflection away from the operative corridor.

Built for reliable anatomical grip, this Dockhorn Retractor configuration features an overall length of 20.3 cm (8 inches). Its working end incorporates a rigid L-shaped blade engineered to comfortably hook beneath fascial borders and muscular planes. The instrument is available in two selectable blade dimensions: 12.0 mm × 40.0 mm (\$64) and 14.0 mm × 55.0 mm (\$69). The specialized concave upper surface stabilizes retracted tissues, preventing slippage under continuous counter-traction.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Dockhorn Retractor to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for outpatient clinics and hospital operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)