



RETRACTOR MAKER

PRODUCT BROCHURE OF

Jolls Friedman Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Jolls Friedman Self Retaining Retractor

SKU RM-JFSRR-1

\$94 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Jolls Friedman Self Retaining Retractor	Dimensions / Details: 6" (15.24 cm) length, fenestrated curved jaws, screw lock	\$94
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging	\$94

PRODUCT DESCRIPTION

The Jolls Friedman Self Retaining Retractor is a specialized endocrine and cervical surgical retractor engineered specifically to provide continuous, hands-free soft-tissue and strap muscle displacement during thyroid and parathyroid procedures. Designed with an anatomical curvature that conforms over the anterior neck, this instrument maintains a wide operative window without requiring manual assistant holding or crowding the surgical field.

Built to clinical framework specifications, this dedicated Jolls Friedman Self Retaining Retractor configuration features an overall instrument length of 6 inches (15.24 cm, \$94). Its working ends incorporate an anatomically curved jaw profile featuring smooth, fenestrated end-surfaces that hold the sternocleidomastoid and strap muscles securely without crushing muscle fibers. The base mechanical locking screw-lock control system allows micro-adjustable tension control that locks firmly under counter-traction without slipping.

Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Jolls Friedman Self Retaining Retractor to withstand repetitive hospital autoclave cycles, providing endocrine, ENT, and general surgical teams worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)