



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hohmann Retractor Rounded End

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hohmann Retractor Rounded End

SKU RM-HRRE-1

\$74 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8 5/8" (21.91 cm) — 19 mm round blade	RM-HRRE-1-1	\$74
9 3/4" (24.76 cm) — 19 mm round blade	RM-HRRE-1-2	\$79
11.0" (27.94 cm) — 19 mm round blade	RM-HRRE-1-3	\$84
11 1/2" (29.21 cm) — 19 mm round blade	RM-HRRE-1-4	\$89
12.0" (30.48 cm) — 19 mm round blade	RM-HRRE-1-5	\$94
12 1/8" (30.80 cm) — 19 mm round blade	RM-HRRE-1-6	\$99

PRODUCT DESCRIPTION

The Hohmann Retractor Rounded End is an essential orthopedic bone leverage instrument engineered specifically to provide safe, atraumatic soft-tissue displacement and joint exposure during major reconstructive procedures. Designed with an ultra-smooth, rounded tip beak, this retractor hooks over long bone cortices cleanly, eliminating the risk of accidental neurovascular puncture or periosteal laceration common with sharp-tipped instruments.

Built for complete anatomical adaptability, the Hohmann Retractor Rounded End is manufactured with a uniform 19 mm round blade width profile across an extensive array of six overall lengths: 8 5/8" (\$74), 9 3/4" (\$79), 11" (\$84), 11 1/2" (\$89), 12" (\$94), and 12 1/8" (\$99). The curved blade architecture cradles bone contours smoothly, while the ergonomic handle provides surgical assistants with sustained, non-slip pulling leverage throughout prolonged joint reaming and femoral preparation.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Hohmann Retractor Rounded End to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for orthopedic joint replacement and trauma centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)