



RETRACTOR MAKER

PRODUCT BROCHURE OF

Richardson Retractor Grip Handle

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Richardson Retractor Grip Handle

SKU RM-RRGH-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
9 ½" (24.13 cm) — Hollow, Lamb, or Lightweight Grip	RM-RRGH-1-1	\$64
10.0" (25.40 cm) — Hollow, Lamb, or Lightweight Grip	RM-RRGH-1-2	\$69
11.0" (27.94 cm) — Hollow, Lamb, or Lightweight Grip	RM-RRGH-1-3	\$74

PRODUCT DESCRIPTION

The Richardson Retractor Grip Handle is an advanced ergonomic refinement of the classic abdominal retractor, engineered specifically to provide deep muscular displacement while reducing user fatigue during major open procedures. Designed with an elongated right-angled concave blade terminating in a soft retaining lip, this instrument allows surgical assistants to apply continuous counter-traction on dense abdominal walls and visceral planes without slipping.

Built to withstand high-frequency autoclave sterilization while remaining completely corrosion resistant, this Richardson Retractor Grip Handle series is highly adaptable across three overall lengths: 9 ½ inches (24.13 cm, \$64), 10 inches (25.40 cm, \$69), and 11 inches (27.94 cm, \$74). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), surgeons can select from three specialized handle profiles: a Hollow Grip Handle profile built with a lightweight framework to minimize strain, a Lamb Handle profile ergonomically molded for steady palm bracing, or a Lightweight Grip Handle profile. Its balanced chassis provides operating assistants with fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Richardson Retractor Grip Handle to withstand heavy mechanical loads and repeated hospital sterilization cycles, providing general, gynecological,

and orthopedic operating rooms worldwide with exceptional clinical durability, corrosion resistance, and surgical reliability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)