



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Perkins Self-Retaining Retractor

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# Perkins Self-Retaining Retractor

**SKU** RM-PSRR-99-1

**\$99 USD**

## AVAILABLE SIZES & PRICING

| SIZE / OPTION                                       | SKU / DETAILS  | PRICE (USD) |
|-----------------------------------------------------|----------------|-------------|
| Standard Solid Prong Pattern — 13 cm, Left or Right | RM-PSRR-99-1-1 | \$99        |
| 2 Prongs & Pointed Blade — 11 cm                    | RM-PSRR-99-1-2 | \$99        |
| 2 Prongs & Pointed Blade — 13 cm                    | RM-PSRR-99-1-3 | \$99        |
| 3 Prongs & Pointed Blade — 11 cm                    | RM-PSRR-99-1-4 | \$99        |
| 3 Prongs & Pointed Blade — 13 cm                    | RM-PSRR-99-1-5 | \$99        |
| Curved Fenestrated Blade — 13 cm, Left or Right     | RM-PSRR-99-1-6 | \$99        |
| Solid Spike Blade — 12.7 cm, Left or Right          | RM-PSRR-99-1-7 | \$99        |

## PRODUCT DESCRIPTION

The Perkins Self-Retaining Retractor is an internationally acclaimed surgical spreader engineered specifically to maintain a stable, hands-free open operative field across vascular, neurological, and orthopedic surgeries. Designed with diverse prong arrays and specialized blade profiles, this instrument optimizes visualization and procedural precision while eliminating manual assistant fatigue during intricate dissections.

Built for complete procedural adaptability, this Perkins Self-Retaining Retractor family (\$99) features ergonomic finger ring handles integrated with a secure mechanical ratchet lock system to ensure steady positioning. The series encompasses five distinct architectural models: a Standard Solid Prong Pattern (13 cm), a 2-Prongs & Pointed Blade Pattern (11 cm or 13 cm), a 3-Prongs & Pointed Blade Pattern (11 cm or 13 cm with security

serrations), a Curved Fenestrated Blade Pattern (13 cm with 3 blunt prongs), and a Solid Spike Blade Pattern (12.7 cm with an 11 × 25 mm footprint, 7 mm spike, and 80 mm maximum spread).

Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Perkins Self-Retaining Retractor to withstand repetitive hospital autoclave cycles, delivering high corrosion resistance, unyielding structural stability, and complete clinical dependability for operating suites worldwide.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)