



RETRACTOR MAKER

PRODUCT BROCHURE OF

Vascular Lobster System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Vascular Lobster System

SKU RM-VLS-1

\$2,699 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Lobster Frame	—
1	Lobster Table Clamp with Vertical Post & Frame	—
8	Lobster Retractor Clamp	—
2	Lobster Mayo Swivel Retractor	$\frac{3}{4}'' \times 2''$
1	Lobster Mayo Swivel Retractor	$3\frac{1}{2}'' \times 2''$
1	Lobster Fence Swivel Retractor	$4'' \times 6''$
1	Lobster Splanchnic Swivel Retractor	$2'' \times 4\frac{1}{2}''$
1	Lobster Splanchnic Swivel Retractor	$2'' \times 6''$
2	Lobster Renal Swivel Retractor	$1'' \times 6''$
1	Lobster Deaver Swivel Retractor	$\frac{3}{4}'' \times 4\frac{1}{2}''$
1	Lobster Malleable Swivel Retractor (perforated)	$2'' \times 6\frac{1}{2}''$
1	Lobster Malleable Swivel Retractor (perforated)	$3'' \times 6\frac{1}{2}''$
1	Lobster Harrington Sweetheart Retractor	$2'' \times 5''$

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Sterilization Box	—

PRODUCT DESCRIPTION

The Vascular Lobster System is an adaptable table-mounted surgical retractor engineered to provide reliable, hands-free exposure during high-risk open vascular interventions. Anchored rigidly to the operating table rail via a vertical field post and mounting clamp, this self-retaining platform eliminates surgical fatigue and prevents accidental blade shifts when operating near fragile major vessels.

Configured specifically for peripheral vascular and abdominal aortic access, the Vascular Lobster System features a low-profile Lobster structural frame and eight fast-locking retractor clamps. The 10-piece swivel blade complement provides targeted holding power, including narrow Renal vein retractors, deep visceral Splanchnic blades (up to 6"), protective Fence blades, Deaver, Harrington Sweetheart, Mayo, and perforated malleable strips that conform to vessel walls without causing shear injury.

Every component is manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a smooth anti-glare satin finish. Retractor Maker supplies the Vascular Lobster System non-sterile and fully autoclavable in a custom-fitted sterilization box, delivering hospital-grade durability and complete traceability for surgical supply chains.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)