



RETRACTOR MAKER

PRODUCT BROCHURE OF

Semb Lung/Muscle Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Semb Lung/Muscle Retractor

SKU RM-SLMR-89-1

\$89 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
25.5 cm Contoured Model — 25.5 cm with Contoured Grip	RM-SLMR-89-1-1	\$89
25.0 cm Flat Handle Model — 25.0 cm with Flat Handle	RM-SLMR-89-1-2	\$89

PRODUCT DESCRIPTION

The Semb Lung/Muscle Retractor is an internationally acclaimed thoracic surgical retractor engineered specifically to provide safe, broad-surface tissue displacement of the lung lobes, chest wall musculature, and deep thoracic structures during open chest procedures. Designed with an anatomical downward-curved blade, this instrument allows surgical assistants to apply sustained counter-traction across heavy tissue planes without causing pulmonary contusions or parenchymal tears.

Built to provide supreme durability and backed by a lifetime warranty against rust and defects, this Semb Lung/Muscle Retractor series offers two primary structural configurations (\$89). The 25.5 cm Length Configuration features a specialized contoured grip handle that gives assistants superior manual leverage and is customizable across three distinct working blade widths: $\frac{1}{2}$ inch wide, $\frac{3}{4}$ inch wide, and 1 inch wide. The alternative 25.0 cm Length Configuration features an integrated flat handle paired with a multi-dimensional blade footprint measuring $\frac{5}{8}$ inches wide $\times$ $1\frac{1}{4}$ inches long. Both designs incorporate rounded, atraumatic blade borders to safeguard vital cardiopulmonary anatomy.

Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Semb Lung/Muscle Retractor to withstand repetitive hospital sterilization cycles, providing cardiothoracic, vascular, and trauma surgical teams worldwide with exceptional structural rigidity, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)