



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cope Retractor Double Ended

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cope Retractor Double Ended

SKU RM-CRDE-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Cope Retractor Double Ended	18 cm overall length, 6.0 mm & 11.0 mm dual blades	\$29
Protective Sleeve / Box	Individual transit packaging	\$29

PRODUCT DESCRIPTION

The Cope Retractor Double Ended is a highly functional handheld surgical instrument engineered specifically to provide atraumatic soft-tissue and wound margin retraction across a broad range of clinical specialties. Designed with contrasting blade dimensions integrated on a single chassis, this retractor allows surgical teams to transition instantly between superficial and intermediate tissue depths without exchanging instruments.

Built for exceptional manual control, this Cope Retractor Double Ended configuration features an overall length of 7 ½ inches (18 cm, \$29). Its working ends incorporate a narrow 6.0 mm wide blade on one end for precise, localized wound retraction and an 11.0 mm wide blade on the opposite end for broader muscle and fascial displacement. The smooth, contoured blade edges ensure reliable tissue traction while preventing localized ischemia or tears. Its balanced central handle provides a fatigue-free grip for assistant staff.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Cope Retractor Double Ended to withstand repetitive hospital sterilization cycles, providing general, plastic, and ENT surgery suites with high corrosion resistance and dependable durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)