



RETRACTOR MAKER

PRODUCT BROCHURE OF

Volkmann Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Volkmann Retractor

SKU RM-VR-SERIES-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
2 Prongs (Blunt / Sharp) — 8 mm - 11 mm	RM-VR-SERIES-1-1	\$64
3 Prongs (Blunt / Sharp) — 10 mm - 18 mm	RM-VR-SERIES-1-2	\$69
4 Prongs (Blunt / Sharp) — 13 mm - 28 mm	RM-VR-SERIES-1-3	\$74
6 Prongs (Blunt / Sharp) — 18 mm - 28 mm	RM-VR-SERIES-1-4	\$79
8 Prongs (Blunt / Sharp) — 28 mm - 42 mm	RM-VR-SERIES-1-5	\$84

PRODUCT DESCRIPTION

The Volkmann Retractor is an internationally acclaimed surgical rake retractor engineered specifically to provide dependable, non-slip skin flap elevation and soft-tissue displacement across general, plastic, and orthopedic surgical procedures. Designed with curved rake prongs, this instrument allows surgical assistants to apply uniform manual counter-traction across cutaneous and subcutaneous planes without causing focal tissue tears or margin ischemia.

Built to withstand high-frequency autoclave sterilization cycles without degrading, this Volkmann Retractor line is exceptionally customizable through a comprehensive matrix of dimensions. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers two overall reaches: an 11.5 cm reach footprint for compact fields and a 22.0 cm reach footprint for deep surgical access. Surgeons can select across five prong configurations: 2 (\$64), 3 (\$69), 4 (\$74), 6 (\$79), or 8 (\$84) prongs, selectable with Blunt prong tips

engineered for safe, extra-atraumatic soft tissue manipulation, or Sharp prong tips engineered to provide an optimal, slip-resistant hold on dense fascia. Blade width options span from 8 mm up to 42 mm wide (including 8, 10, 11, 13, 18, 19, 28, and 42 mm).

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Volkmann Retractor to uphold strict clinical standards, providing operating theaters, hand trauma suites, and reconstructive clinics worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)