



RETRACTOR MAKER

PRODUCT BROCHURE OF

Doyen Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Doyen Vaginal Retractor

SKU RM-DVR-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Doyen Vaginal Retractor	Handheld broad-blade vaginal retractor with ring handle	\$79
Protective Sleeve / Box	Individual transit packaging	\$79

PRODUCT DESCRIPTION

The Doyen Vaginal Retractor is a foundational gynecological surgical instrument engineered specifically to provide broad, atraumatic displacement of the vaginal walls during pelvic and vaginal operative procedures. Designed with an anatomically curved distal paddle, this retractor depresses the posterior or anterior vaginal mucosa smoothly, providing the operating surgeon with an unhindered visual corridor to the cervix and fornices.

Built for superior pelvic exposure, this Doyen Vaginal Retractor configuration (\$79) incorporates a broad blade face terminating in a rounded retaining lip that prevents mucosal slippage without causing localized ischemic injury. Its ergonomically angled handle features a finger ring grip that gives surgical assistants exceptional mechanical advantage, ensuring stable counter-traction throughout lengthy vaginal reconstructions.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Doyen Vaginal Retractor to withstand repetitive hospital sterilization cycles, providing gynecological, obstetrical, and urogynecological surgery suites worldwide with exceptional durability, corrosion resistance, and clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)