



RETRACTOR MAKER

PRODUCT BROCHURE OF

Deep Gelpi Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Deep Gelpi Self Retaining Retractor

SKU RM-DGSRR-1

\$129 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
7 ½" (approx. 19.0 cm) — Angled chassis, sharp out-turned prongs — Standard lumbar spine laminectomy and joint exposure	RM-DGSRR-1-1	\$129
10 ½" (approx. 26.7 cm) — Extended angled chassis, sharp out-turned prongs — Deep pelvic, bariatric, and multi-level spinal fusion	RM-DGSRR-1-2	\$139

PRODUCT DESCRIPTION

The Deep Gelpi Self Retaining Retractor is a specialized surgical retractor engineered specifically to provide powerful, hands-free soft-tissue and muscular retraction during deep spinal, orthopedic, and pelvic operations. Designed with out-turned pointed prongs on an angled chassis, this retractor embeds securely into dense fascial edges and paraspinal muscle, establishing an unyielding operative corridor without manual assistance.

Built for deep-cavity exposure, this Deep Gelpi Self Retaining Retractor configuration features an angled working chassis that keeps the finger-ring handles flat against patient contours, preventing interference with surgical navigation tools and microscopes. Surgeons can select between two working lengths: 7 ½ inches (\$129) for standard spinal corridors or an extended 10 ½ inches (\$139) reach for deep pelvic and bariatric lumbar procedures. Its heavy-duty mechanical ratchet lock maintains continuous tension without slipping.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Deep Gelpi Self Retaining Retractor to withstand heavy muscular counter-traction, providing spine surgeons, neurosurgeons, and orthopedic trauma teams with exceptional durability and complete autoclave compatibility.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)