



RETRACTOR MAKER

PRODUCT BROCHURE OF

Tebbetts Double-Ended Breast Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Tebbetts Double-Ended Breast Retractor

SKU RM-TDEBR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
23 cm Reach Layout — 37×70 mm long blade / 22×54 mm short blade — Standard inframammary augmentation and periareolar pocket dissection	RM-TDEBR-1-1	\$69
31 cm (12.25") Reach — 37×70 mm long blade / 22×54 mm short blade — Deep subpectoral pocket creation, transaxillary augmentation, and revision mammoplasty	RM-TDEBR-1-2	\$89

PRODUCT DESCRIPTION

The Tebbetts Double-Ended Breast Retractor is a specialized double-ended aesthetic surgical retractor engineered specifically to provide versatile, multi-depth tissue elevation and pocket exposure during breast augmentation and reconstruction. Designed with two distinct curved blades on opposing ends, this instrument allows plastic surgeons to transition effortlessly between superficial incision retraction and deep subpectoral pocket exploration without exchanging tools.

Built for superior procedural agility, this Tebbetts Double-Ended Breast Retractor features a specialized gold-plated handle to maximize maneuverability and tactile grip control during lengthy operative cases. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument is available across two overall longitudinal reaches: a 23 cm reach layout (\$69) for standard access and an extended 31 cm (12.25 inches) reach layout (\$89) for deep pocket visualization. The dual-ended blade footprint dimensions incorporate a Longer Blade measuring 37 mm wide × 70 mm long on one end paired with a Shorter Blade measuring 22 mm wide × 54 mm long on the reverse end, both contoured with smooth atraumatic margins.

Forged with an anti-glare satin matte finish on working blades and an elegant gold-plated handle, Retractor Maker manufactures the Tebbetts Double-Ended Breast Retractor to withstand repetitive hospital autoclave cycles, delivering exceptional corrosion resistance, structural rigidity, and complete clinical durability for aesthetic and reconstructive plastic surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)