



RETRACTOR MAKER

PRODUCT BROCHURE OF

Omni Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Omni Retractor System

SKU RM-OR-1

\$2,499 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Omni Clamp Sterile Field Post	Heavy-Duty table rail clamp and post
1	Standard Wishbone Frame	Split fork-style structural ring frame
1	Hinged Extension Arm	19" (48.3 cm) for extended reach and flexibility
8	SnapClamp Blade Holders	Quick-connect, slide-on frame clamps
2	Mayo Swivel Retractors	2¾" × 2" (7.0 × 5.1 cm)
1	Mayo Swivel Retractor	3½" × 2" (8.9 × 5.1 cm)
1	Fence Swivel Retractor	4" × 6" (10.2 × 15.2 cm) broad muscle wall blade
1	Splanchnic Swivel Retractor	2" × 4½" (5.1 × 11.4 cm) deep visceral blade
1	Splanchnic Swivel Retractor	2" × 6" (5.1 × 15.2 cm) deep visceral blade
2	Renal Vein Swivel Retractors	1" × 6" (2.5 × 15.2 cm) narrow vascular blades
1	Deaver Swivel Retractor	2¾" × 4½" (7.0 × 11.4 cm) curved blade
1	Malleable Slotted Swivel Retractor	2" × 6½" (5.1 × 16.5 cm) bendable strip
1	Malleable Slotted Swivel Retractor	3" × 6½" (7.6 × 16.5 cm) bendable strip
1	Harrington (Sweetheart) Swivel Retractor	2" × 5" (5.1 × 12.7 cm) heart-shaped blade

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Storage Box	Wooden storage case with customized padding slots

PRODUCT DESCRIPTION

The Omni Retractor System is an advanced, table-mounted self-retaining retractor engineered to deliver customizable, low-profile operative exposure during complex open surgical procedures. Designed around an adaptable Wishbone split-frame, a heavy-duty field post clamp, a 19-inch hinged extension arm, and eight quick-connect SnapClamp blade holders, this system delivers dependable multi-directional traction without occupying vital space around the sterile field.

Engineered for precision across vascular, transplant, general, urological, and deep pelvic procedures, the Omni Retractor System features an extensive complement of specialized swivel blades—including Mayo, Fence, Splanchnic, narrow Renal Vein, Deaver, Harrington (Sweetheart), and malleable slotted retractors. Every component is machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin finish, and the full assembly is housed in a custom wooden storage box for safe transit and organized facility storage.

Manufactured in our dedicated production facility, Retractor Maker supplies the Omni Retractor System directly to hospitals, surgical centers, and healthcare distributors worldwide, backed by complete material traceability and direct factory support.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)