



RETRACTOR MAKER

PRODUCT BROCHURE OF

Lateral Nonconductive Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Lateral Nonconductive Vaginal Retractor

SKU RM-LNCVR-63-1

\$69 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Open Shanks Profile	Blade Width: 6.3 cm · Primary Clinical Role: Enhanced intraoperative instrument clearance during LEEP/LLETZ · Price (USD): \$69.00
1	Closed Elongated Shanks Profile	Blade Width: 6.3 cm · Primary Clinical Role: Streamlined lateral vaginal mucosal displacement and laser ablation · Price (USD): \$69.00
1	Protective Sleeve / Box	Blade Width: Standard Transit Packaging · Primary Clinical Role: — · Price (USD): Included

PRODUCT DESCRIPTION

The Lateral Nonconductive Vaginal Retractor is a specialized gynecological surgical instrument engineered specifically to provide safe lateral vaginal wall retraction during loop electrosurgical excision procedures (LEEP) and laser ablation surgery. Coated with a high-integrity, non-conductive dielectric polymer layer, this instrument prevents stray electrical currents and thermal transfer from damaging delicate collateral vaginal tissues.

Built in a standard configuration layout, this Lateral Nonconductive Vaginal Retractor model incorporates a 6.3 cm wide blade footprint (\$69). Surgeons can customize the instrument by selecting between Open shanks for expanded instrument maneuverability or Closed elongated shanks profiles for streamlined insertion. The surface profile combines highly polished handling interfaces with a specialized insulated working layer that withstands high-frequency electrosurgical cutting and coagulation without degradation.

Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin finish on exposed metal elements, Retractor Maker produces the Lateral Nonconductive Vaginal Retractor to strict clinical tolerances, delivering exceptional procedural hygiene, high dielectric resistance, and complete

autoclave durability for hospital operating theaters and women's health clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)