



RETRACTOR MAKER

PRODUCT BROCHURE OF

Kocher Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Kocher Retractor

SKU RM-KCR-1

\$49 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	20.3 cm layout with 20 mm × 12 mm blade	Primary Clinical Role: Superficial wound margin and thyroid lobe elevation · Price (USD): \$49.00
1	21 cm layout with 8×25, 12×40, 12×45 mm blades	Primary Clinical Role: Intermediate soft-tissue retraction and appendectomy · Price (USD): \$54.00
1	21 cm layout with 12×60, 12×65 mm blades	Primary Clinical Role: Deep fascial retraction and extremity muscle holding · Price (USD): \$54.00
1	23 cm layout with 18×40, 20×38, 20×60 mm blades	Primary Clinical Role: Standard adult laparotomy and orthopedic joint exposure · Price (USD): \$64.00
1	23 cm layout with 25×60, 30×75, 40×75 mm blades	Primary Clinical Role: Broad muscular displacement in pelvic and hip surgery · Price (USD): \$64.00
1	Protective Sleeve / Box	Primary Clinical Role: Standard Transit Packaging · Price (USD): Included

PRODUCT DESCRIPTION

The Kocher Retractor is an essential handheld surgical retractor engineered specifically to provide dependable, deep soft-tissue displacement and muscle holding during general, thyroid, and orthopedic procedures. Designed with an elongated right-angled blade terminating in a soft retaining lip, this instrument allows surgical assistants to maintain constant, non-slip counter-traction on thick tissue layers without crushing underlying structures.

Built for complete procedural adaptability, this Kocher Retractor series maintains a fixed finger-ring layout while featuring highly adaptive size menus combining instrument length and blade surface footprint dimensions. The

instrument is available in four overall lengths: 20.3 cm (\$49), 21 cm (\$54), 22.8 cm (\$59), and 23 cm (\$64). Surgeons can select from 12 distinct combined configurations, including 20.3 cm with a 20 × 12 mm blade, 21 cm layouts with blades from 8 × 25 mm to 12 × 65 mm, and 23 cm layouts with broad blades spanning 18 × 40 mm up to 40 × 75 mm.

Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Kocher Retractor to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for operating rooms worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)