



RETRACTOR MAKER

PRODUCT BROCHURE OF

Henly Retractor Frame

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Henly Retractor Frame

SKU RM-HRF-6CM-1

\$124 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Henly Retractor Frame	6 cm spread capacity, finger-ring ratchet chassis, center blade socket	\$124
Protective Sleeve / Box	Individual transit packaging	\$124

PRODUCT DESCRIPTION

The Henly Retractor Frame is a foundational self-retaining surgical retractor chassis engineered specifically to provide reliable, multi-directional soft-tissue exposure during vascular, general, and orthopedic procedures. Built to host modular side and center blade attachments, this mechanical framework allows surgeons to convert standard two-arm lateral spreading into comprehensive three-way operative access without requiring assistant hands.

Built to clinical framework standards, this Henly Retractor Frame configuration features a maximum spread capacity of 6 cm (60 mm total arm separation width, \$124) that accommodates variable wound dimensions. Its handle assembly incorporates a dual finger ring layout with built-in multi-position locking ratchet gears that sustain continuous retraction tension automatically. The integrated center crossbar accommodates interchangeable center blades, ensuring rigid stability throughout delicate vascular dissections.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Henly Retractor Frame to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, smooth mechanical gear tracking, and complete autoclave durability for operating rooms worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)