



RETRACTOR MAKER

PRODUCT BROCHURE OF

Converse Alar Retractor Double Ended

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Converse Alar Retractor Double Ended

SKU RM-CARDE-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
10 mm Wide Opening	RM-CARDE-1-1	\$49
14 mm Wide Opening	RM-CARDE-1-2	\$54

PRODUCT DESCRIPTION

The Converse Alar Retractor Double Ended is a classic rhinoplasty surgical retractor engineered specifically to provide safe, non-traumatic retraction of the nostril margin and alar cartilage during aesthetic and reconstructive nasal procedures. Designed with an ergonomic double-ended profile, this instrument allows surgeons to adapt instantly to variable nostril diameters without changing tools.

Built for precision handling in tight surgical spaces, this Converse Alar Retractor Double Ended model features a compact 4-inch (10.2 cm) total profile. Its working ends incorporate specialized curved hook blades featuring selectable dual-ended openings: a narrow 10 mm wide opening on one end (\$49) for tight vestibular corridors and an expanded 14 mm wide opening on the opposite end (\$54) for broader alar base elevation. The smooth, rounded blade borders prevent tearing of delicate vestibular skin and mucosal linings.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Converse Alar Retractor Double Ended to strict clinical standards, delivering high corrosion resistance, excellent tactile sensitivity, and complete autoclave durability for plastic and ENT operating suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)