



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Coryllos Retractors

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



## Coryllos Retractors

**SKU** RM-COR-THOR-1

**\$99 USD**

### AVAILABLE SIZES & PRICING

| SIZE / OPTION               | SKU / DETAILS   | PRICE (USD) |
|-----------------------------|-----------------|-------------|
| 22 mm × 82 mm — 9 ½" - 10"  | RM-COR-THOR-1-1 | \$99        |
| 43 mm × 108 mm — 9 ½" - 10" | RM-COR-THOR-1-2 | \$109       |
| 65 mm × 115 mm — 9 ½" - 10" | RM-COR-THOR-1-3 | \$119       |

### PRODUCT DESCRIPTION

Coryllos Retractors are specialized cardiothoracic surgical instruments engineered specifically to provide deep, uncompromised exposure of the thoracic cavity and chest wall during thoracotomy and lung resection procedures. Designed with an expansive, anatomically curved blade, these instruments allow surgical assistants to depress lung parenchyma and retract intercostal musculature smoothly without causing focal tissue contusions.

Built for heavy thoracic leverage, the Coryllos Retractors series features an operational reach spanning 9 ½ inches to 10 inches (approx. 24–25.4 cm) to provide excellent physical leverage. The line is manufactured in three standardized blade dimensions: 22 mm × 82 mm (\$99), 43 mm × 108 mm (\$109), and an extra-wide 65 mm × 115 mm (\$119) profile. Its teardrop-aperture ring handle gives surgical assistants a secure, non-slip grip that significantly reduces hand strain during multi-hour open thoracic procedures.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures Coryllos Retractors to withstand heavy mechanical loads and repeated hospital autoclaving cycles, providing cardiothoracic and general surgical teams with long-lasting clinical dependability.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)