



RETRACTOR MAKER

PRODUCT BROCHURE OF

Heaney Hysterectomy Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Heaney Hysterectomy Retractor

SKU RM-HEHR-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
44 mm × 25 mm — 17 cm or 27 cm	RM-HEHR-1-1	\$84
70 mm × 25 mm — 17 cm or 27 cm	RM-HEHR-1-2	\$84
73 mm × 28.5 mm — 17 cm or 27 cm	RM-HEHR-1-3	\$89
114 mm × 25 mm — 17 cm or 27 cm	RM-HEHR-1-4	\$89

PRODUCT DESCRIPTION

The Heaney Hysterectomy Retractor is a foundational gynecological surgical instrument engineered specifically to provide broad, non-slip retraction of the vaginal wall and pelvic tissues during vaginal hysterectomy and pelvic reconstruction. Designed with an elongated right-angled blade terminating in a pronounced retaining lip, this instrument allows surgeons to maintain deep operative visualization of the cervix, uterosacral ligaments, and vaginal fornices.

Built for complete procedural adaptability, this Heaney Hysterectomy Retractor series is available in two overall lengths: a compact 17 cm model (\$84) for shallow fields and an extended 27 cm model (\$89) for deep pelvic access. Surgeons can select across four distinct blade profiles: 44 mm × 25 mm, 70 mm × 25 mm, 73 mm × 28.5 mm, and an extra-deep 114 mm × 25 mm configuration. Its ergonomic handle delivers steady manual leverage, preventing tissue slippage without lacerating delicate vaginal mucosa.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Heaney Hysterectomy Retractor to withstand repetitive hospital sterilization cycles, providing obstetrical, gynecological, and urogynecological teams worldwide with exceptional durability, corrosion resistance, and complete autoclave compatibility.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)