



RETRACTOR MAKER

PRODUCT BROCHURE OF

St. Marks Pelvic Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



St. Marks Pelvic Retractor

SKU RM-SMPR-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
11 ½" Lipped Variant — 11 ½" (29.2 cm)	RM-SMPR-1-1	\$79
13" Lipped Variant — 13" (33.0 cm)	RM-SMPR-1-2	\$84
13" Non-Lipped Variant — 13" (33.0 cm)	RM-SMPR-1-3	\$89

PRODUCT DESCRIPTION

The St. Marks Pelvic Retractor is an internationally acclaimed deep-cavity surgical instrument engineered specifically to provide safe, expansive visceral displacement and muscular retraction during deep colorectal, pelvic, and gynecological procedures. Designed with an elongated right-angled blade, this instrument allows surgical teams to displace the bladder, uterus, and pelvic peritoneum cleanly, maintaining an unobstructed visual corridor into the presacral space and deep pelvic floor.

Built to withstand high-frequency autoclave sterilization cycles while providing a lifetime warranty against rust and material defects, this St. Marks Pelvic Retractor series offers three distinct configuration models: an 11 ½-inch Length Variant (\$79) equipped with a Lipped blade measuring 35 mm × 45 mm, a 13-inch Lipped Length Variant (\$84) featuring a 60 mm × 45 mm blade, and a 13-inch Non-Lipped Length Variant (\$89) with a 60 mm × 45 mm blade. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the rounded blade margins prevent visceral perforation, while the robust ergonomic handle delivers fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the St. Marks Pelvic Retractor to withstand immense mechanical counter-traction, providing colorectal surgical suites, urological operating theaters, and gynecologic oncology centers worldwide with unyielding structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)