



RETRACTOR MAKER

PRODUCT BROCHURE OF

Collin Hartmann S-Shaped Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Collin Hartmann S-Shaped Retractor

SKU RM-CHSR-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Collin Hartmann S-Shaped Retractor	16 cm overall length, 52×33 mm & 30×21 mm dual ends	\$29
Protective Sleeve / Box	Individual transit packaging	\$29

PRODUCT DESCRIPTION

The Collin Hartmann S-Shaped Retractor is an adaptable double-ended surgical retractor engineered to provide atraumatic soft-tissue displacement across a wide spectrum of general, gynecological, and reconstructive procedures. Designed with an ergonomic S-shaped curvature, this instrument allows surgical assistants to retract tissue layers from multiple angles without obstructing the primary surgeon's working field.

Built for dual-purpose procedural utility, this Collin Hartmann S-Shaped Retractor model features an overall length of 16 cm (\$29) integrating two distinct working ends on a single chassis. The Large Working Surface End measures 52 mm × 33 mm, optimized for retracting broad abdominal walls, adipose tissue, and muscular flaps. The opposite Small Working Surface End measures 30 mm × 21 mm, engineered for smaller incisions, pediatric fields, and fine anatomical contours.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the Collin Hartmann S-Shaped Retractor non-sterile and fully autoclavable inside protective packaging, delivering high corrosion resistance, unyielding structural rigidity, and complete clinical dependability for operating rooms worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)