



RETRACTOR MAKER

PRODUCT BROCHURE OF

Barsky Nasal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Barsky Nasal Retractor

SKU RM-BNR-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Barsky Nasal Retractor	3 ¾" (9.5 cm) length, 10 × 55 mm blade	\$64
Protective Sleeve / Box	Individual transit packaging	\$64

PRODUCT DESCRIPTION

The Barsky Nasal Retractor is a specialized handheld plastic and reconstructive surgery retractor engineered to provide gentle, broad-surface tissue retraction during rhinoplasty and cleft palate repairs. Designed with an elongated blade profile, this instrument distributes holding pressure evenly across the nasal wall, preventing localized skin ischemia and mucosal tearing during delicate dissections.

Built for delicate facial procedures, this Barsky Nasal Retractor configuration features an overall length of 3 ¾ inches (9.5 cm) with a 10 mm wide by 55 mm long blade. Its broad surface area allows surgical teams to securely elevate the lateral nasal cartilage and alar base, establishing a clear view of the lower lateral cartilages and vestibular vaults without obstructing surgical instruments.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Barsky Nasal Retractor to strict clinical standards, providing plastic and ENT surgical suites with a fully autoclavable, corrosion-resistant instrument ready for daily operative use.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)