



RETRACTOR MAKER

PRODUCT BROCHURE OF

Richardson Retractor Loop Handle

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Richardson Retractor Loop Handle

SKU RM-RRLH-1

\$54 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
9 ½" (24.13 cm) — Right-angle curved blade, loop handle	RM-RRLH-1-1	\$54
10.0" (25.40 cm) — Right-angle curved blade, loop handle	RM-RRLH-1-2	\$59

PRODUCT DESCRIPTION

The Richardson Retractor Loop Handle is an essential handheld surgical retractor engineered specifically to provide reliable soft-tissue and muscular displacement during open abdominal and pelvic operations. Designed with a classic right-angled curved blade terminating in a downward retaining lip, this instrument allows surgical assistants to apply steady manual counter-traction across heavy tissue planes without causing tissue slippage or local margin ischemia.

Built to withstand heavy physical demands and frequent autoclave sterilization while remaining completely corrosion resistant, this Richardson Retractor Loop Handle series is offered in two overall lengths: 9 ½ inches (24.13 cm, \$54) and 10 inches (25.40 cm, \$59). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry layout incorporates a single-ended right-angle curved blade that holds muscle bellies and viscera securely. The handle assembly features an integrated ergonomic Loop Handle framework tailored to minimize manual hand fatigue and guarantee a non-slip grip during multi-hour surgical procedures.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Richardson Retractor Loop Handle to withstand repetitive hospital autoclave cycles, providing general, gynecological, and trauma surgical operating theaters worldwide with unyielding structural rigidity, high corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)