



RETRACTOR MAKER

PRODUCT BROCHURE OF

Lone Star Retractor (General, Head & Neck Surgeries)

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Lone Star Retractor (General, Head & Neck Surgeries)

SKU RM-LSR-GHN-1

\$499 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small Ring — 22.0 cm x 11.3 cm (General/Head & Neck)	RM-LSR-GHN-1-1	\$484
Standard Ring — 28.3 cm x 18.3 cm (General/Head & Neck)	RM-LSR-GHN-1-2	\$499
Large Ring — 32.5 cm x 18.3 cm (General/Head & Neck)	RM-LSR-GHN-1-3	\$514

PRODUCT DESCRIPTION

The Lone Star Retractor (General, Head & Neck Surgeries) is a highly adaptable, self-retaining retractor frame engineered to provide clear, hands-free visualization during delicate soft-tissue dissections. By holding surgical wound margins open with multi-directional elastic stay hooks anchored into its perimeter notches, this system eliminates manual retraction fatigue and minimizes operative crowding.

Designed specifically for anatomical adaptability in tight surgical corridors, this Lone Star Retractor (General, Head & Neck Surgeries) configuration conforms smoothly to uneven bodily contours around the neck, axilla, and abdominal wall. The frame's precision-milled slot geometry allows surgical teams to position and tension elastic stays incrementally, providing uniform exposure and minimizing focal tissue ischemia along wound margins.

Manufactured from high-strength anodized aluminum and Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), Retractor Maker supplies the Lone Star Retractor (General, Head & Neck Surgeries) non-sterile and fully autoclavable inside an individual protective box, giving surgical centers long-term reliability and effortless intraoperative adjustability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)