



RETRACTOR MAKER

PRODUCT BROCHURE OF

Swivel Ring Nasal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Swivel Ring Nasal Retractor

SKU RM-SRNR-11-1

\$59 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Swivel Ring Nasal Retractor	11 cm overall length, 8 mm wide blade, 360° swivel ring	\$59
Protective Sleeve / Box	Individual transit packaging	\$59

PRODUCT DESCRIPTION

The Swivel Ring Nasal Retractor is an innovative handheld rhinoplasty and facial plastic surgical retractor engineered specifically to provide dynamic, multi-directional elevation of the nasal skin envelope and lower lateral cartilages. Designed with a 360-degree rotating finger ring, this instrument allows operating assistants to adjust the manual pulling angle continuously without releasing tissue contact or occluding the surgeon's operative sightline.

Built to ensure a clean, contained operating field, this Swivel Ring Nasal Retractor configuration features an overall instrument length of 11 cm paired with an 8 mm wide blade (\$59). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to ensure lasting corrosion resistance, its working end incorporates a smoothly contoured curved blade that cradles delicate alar cartilage safely without causing localized crush injury or mucosal tears. The swivel ring mechanism rotates smoothly under tension, providing ergonomic manual control during intricate dorsal remodeling.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Swivel Ring Nasal Retractor to withstand repetitive hospital autoclave cycles, providing cosmetic, ENT, and facial reconstructive surgery suites worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)