



RETRACTOR MAKER

PRODUCT BROCHURE OF

Chandler Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Chandler Retractor

SKU RM-CHR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small (16.0 mm leaf)	RM-CHR-1-1	\$69
Medium (20.0 mm leaf)	RM-CHR-1-2	\$74
Large (25.0 mm leaf)	RM-CHR-1-3	\$79

PRODUCT DESCRIPTION

The Chandler Retractor is an indispensable orthopedic surgical instrument engineered specifically to provide reliable bone leverage and soft-tissue displacement during total hip arthroplasty and reconstructive joint surgery. Designed with an arched leaf-shaped blade that curves backward, this instrument wraps smoothly around bony contours, allowing surgical assistants to maintain continuous counter-traction without slipping off cortical borders.

Built for heavy lower-extremity leverage, the Chandler Retractor series features an overall length of 9.25 inches (23.49 cm) with three selectable tip sizes spanning a broad functional width range: Small (\$69), Medium (\$74), and Large (\$79), with leaf widths measuring from 16.0 mm to 25.0 mm. The blade's rounded blunt edge safeguards delicate femoral vessels, nerves, and periosteum from sharp crush injuries while providing an expansive visual corridor into deep articular spaces.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Chandler Retractor to endure heavy orthopedic leverage without bending, providing joint replacement and trauma teams worldwide with exceptional clinical durability and autoclave longevity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)