



RETRACTOR MAKER

PRODUCT BROCHURE OF

Caspar Lumbar Microdiscectomy Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Caspar Lumbar Microdiscectomy Retractor System

SKU RM-CLMRS-1

\$799 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Lateral Blade	40 mm
1	Lateral Blade	45 mm
1	Lateral Blade	50 mm
1	Lateral Blade	55 mm
1	Lateral Blade	60 mm
1	Lateral Blade	65 mm
1	Lateral Blade	75 mm
1	Specula	40 mm
1	Specula	45 mm
1	Specula	50 mm
1	Specula	55 mm
1	Specula	60 mm
1	Specula	65 mm
1	Specula	75 mm

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Blade Ejector Forceps	Quick-release blade exchange forceps
1	Caspar Microdiscectomy Retractor Frame	Counter-retraction frame
1	Sterilization Box	Aluminium container with organizing tray

PRODUCT DESCRIPTION

The Caspar Lumbar Microdiscectomy Retractor System is a classic, highly regarded surgical retractor platform engineered specifically to provide unobstructed micro-exposure during lumbar disc and canal decompressions. By combining a counter-retraction frame with interchangeable specula and lateral blades, this system establishes a narrow, stable corridor that frees the operative corridor from bulky instrumentation.

Designed for microdiscectomy and focused laminectomies, the Caspar Lumbar Microdiscectomy Retractor System includes 17 precision functional components. The set incorporates a Caspar microdiscectomy retractor frame, seven lateral blades ranging from 40 mm to 75 mm depths, seven matching specula ranging from 40 mm to 75 mm, and dedicated blade ejector forceps for fast, secure intraoperative blade swapping without tissue trauma.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin finish, Retractor Maker supplies the complete Caspar Lumbar Microdiscectomy Retractor System non-sterile and fully autoclavable inside a heavy-duty aluminum sterilization container with dedicated organizing racks, ensuring rapid turnover and long-term durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)