



RETRACTOR MAKER

PRODUCT BROCHURE OF

Talei Necklift Retractor "Benetrator 2"

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Talei Necklift Retractor "Benetrator 2"

SKU RM-TNRB2-9IN-1

\$159 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Talei Necklift Retractor "Benetrator 2"	9" (22.86 cm) overall length, broad curved blade, integrated suction tube	\$159
Protective Sleeve / Box	Individual transit packaging	\$159

PRODUCT DESCRIPTION

The Talei Necklift Retractor "Benetrator 2" is an advanced aesthetic surgical instrument engineered specifically to provide expansive soft-tissue elevation and continuous fluid management during deep plane necklift, platysmaplasty, and submentoplasty procedures. Designed with an anatomical broad blade, this instrument allows plastic surgeons to undermine and elevate thick cervical cutaneous flaps cleanly away from the platysma muscle without compressing the marginal mandibular or cervical branches of the facial nerve.

Built to withstand frequent autoclave sterilization cycles while providing a lifetime warranty against rust, this Talei Necklift Retractor "Benetrator 2" configuration features an overall tool length of 9 inches (22.86 cm, \$159). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry layout incorporates a broad, curved blade featuring an atraumatic flat edge that cradles delicate vascularized flaps securely without causing focal ischemic necrosis. The single-ended design is integrated into an angled grip framework paired with an integrated built-in suction tube channel that clears blood and electrosurgical plume instantly from deep submental tunnels.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Talei Necklift Retractor "Benetrator 2" to uphold the highest standards of facial reconstructive surgery, providing cosmetic operating rooms and outpatient surgery clinics worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)