



RETRACTOR MAKER

PRODUCT BROCHURE OF

Modified Fukuda Type Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Modified Fukuda Type Retractors

SKU RM-MFTR-SHLD-1

\$89 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Modified Fukuda Type Retractor	Handheld orthopedic humeral leverage retractor, fenestrated ring blade	\$89
Protective Sleeve / Box	Individual transit packaging	\$89

PRODUCT DESCRIPTION

Modified Fukuda Type Retractors are indispensable orthopedic surgical instruments engineered specifically to provide dependable, controlled displacement of the humeral head away from the glenoid rim during shoulder reconstruction. Designed with an anatomically curved ring-fenestrated blade, these retractors cradle the proximal humerus securely, allowing surgeons to prepare the glenoid surface without risk of instrument slippage or cartilage shearing.

Built for heavy joint leverage, this Modified Fukuda Type Retractors model (\$89) incorporates an optimized shank curvature that drops the assistant's hands safely below the operative plane. The fenestrated blade architecture allows clear visualization of underlying articular margins while distributing counter-traction forces evenly across the humeral neck cortex. Its ergonomic handle delivers non-slip manual control and maximum leverage during joint manipulation.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures Modified Fukuda Type Retractors to endure heavy orthopaedic torque and repeated hospital autoclave cycles, providing shoulder reconstruction and sports medicine teams worldwide with unyielding clinical durability and performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)