



RETRACTOR MAKER

PRODUCT BROCHURE OF

Aufricht Walter Nasal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Aufricht Walter Nasal Retractor

SKU RM-AWNR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Aufricht Walter Nasal Retractor	13.5 cm overall length, 42 × 7 mm double-ended blades	\$69
Protective Sleeve / Box	Individual transit packaging	\$69

PRODUCT DESCRIPTION

The Aufricht Walter Nasal Retractor is a dual-ended handheld rhinoplasty retractor engineered to provide exceptional surgical versatility and precise soft-tissue traction during cosmetic and reconstructive nasal procedures. By integrating two specialized working ends on a single balanced handle, this instrument allows surgeons to adapt to variable nasal vault depths without changing instruments.

Designed for refined procedural control, this Aufricht Walter Nasal Retractor model features an overall length of 13.5 cm with working blades measuring 42 mm long by 7 mm wide. The double-ended blade architecture provides differing blade angles and curvature profiles, allowing surgical teams to gently support the alar cartilage with one end or elevate the dorsal skin-soft tissue envelope with the opposing end without risking mucosal tears.

Manufactured from high-grade Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker produces the Aufricht Walter Nasal Retractor to withstand rigorous hospital autoclaving cycles, providing plastic and ENT surgical suites worldwide with high durability and surgical precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)