



RETRACTOR MAKER

PRODUCT BROCHURE OF

Black Belt Cervical Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Black Belt Cervical Retractor System

SKU RM-BBCRS-1

\$2,499 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Retractor Frame (Right)	Right-side frame configuration
1	Retractor Frame (Left)	Left-side frame configuration
1	Retractor Frame (Offset)	Offset frame configuration
2	Blade Insertion Clamps	Components used for blade insertion and attachment
14	Pronged Blades	34 mm – 70 mm pronged blade options in multiple sizes
14	Blunt Blades	34 mm – 70 mm blunt blade options in multiple sizes
12	Two-Level Pronged Blades	40 mm – 55 mm two-level pronged blade options
6	Two-Level Blunt Blades	40 mm – 55 mm two-level blunt blade options
1	Sterilization Container	Container for organization, storage, and reprocessing

PRODUCT DESCRIPTION

The Black Belt Cervical Retractor System is an advanced spinal surgical retractor platform engineered to provide secure, stable soft-tissue displacement during anterior cervical spine surgeries. Built to withstand continuous counter-traction without shifting, this system provides unmatched visualization of vertebral bodies and intervertebral disc spaces while protecting vulnerable anatomical structures including the carotid sheath, trachea, and esophagus.

Engineered for single and multi-level spinal interventions, this Black Belt Cervical Retractor System configuration includes three distinct retractor frames (right, left, and offset) and two precision blade insertion clamps. The set comes furnished with an extensive 46-blade inventory comprising 14 pronged blades (34–70 mm), 14 blunt tissue blades (34–70 mm), 12 two-level pronged blades (40–55 mm), and 6 two-level blunt blades (40–55 mm) to provide anatomical versatility across diverse patient physiques.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a non-reflective satin finish, Retractor Maker supplies the Black Belt Cervical Retractor System non-sterile and fully autoclavable in a custom heavy-duty sterilization container with dedicated organizing brackets, ensuring rapid turnover and maximum procedural confidence.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)