



RETRACTOR MAKER

PRODUCT BROCHURE OF

Iron Intern Hydra Double Arm

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Iron Intern Hydra Double Arm

SKU RM-IIHDA-SERIES-1

\$1,699 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
22" Standard Hydra Model	RM-IIHDA-1-1	\$1,699
22" Jumbo Heavy Duty	RM-IIHDA-2-1	\$1,799
26" Jumbo Heavy Duty	RM-IIHDA-3-1	\$1,999

PRODUCT DESCRIPTION

The Iron Intern Hydra Double Arm is an advanced dual-articulated table-mounted surgical retractor system engineered specifically to provide comprehensive bilateral visceral retraction and wide multi-quadrant abdominal wall exposure during major general, bariatric, and transplant surgery. Designed with two fully independent articulated arms radiating from a central table-mounted yoke, this system allows surgical teams to retract opposing wound margins and separate internal organs simultaneously, eliminating the need for multiple manual assistants.

Built to clinical framework specifications, this Iron Intern Hydra Double Arm series features high-load operational configurations across two core reach lengths. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), options include a 22-inch (56 cm) Standard configuration (\$1,699) for adult open laparoscopy and routine general surgery; a 22-inch Jumbo Heavy Duty configuration (\$1,799) for high-load bariatric procedures and dense pelvic layouts; and an extended 26-inch (66 cm) Jumbo Heavy Duty configuration (\$1,999) engineered for complex multi-quadrant transplants, extensive oncology, and extreme exposure conditions. Each arm incorporates independent multi-directional locking mechanics to secure retractor blades rigidly in any anatomical vector.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Iron Intern Hydra Double Arm to withstand heavy surgical counter-traction and repetitive hospital autoclave cycles, providing transplant

operating theaters, bariatric surgical centers, and trauma suites worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)