



RETRACTOR MAKER

PRODUCT BROCHURE OF

Mayo Abdominal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Mayo Abdominal Retractor

SKU RM-MAR-28-1

\$109 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard — 11 in, 70 mm curved blade	RM-MAR-28-1-1	\$109

PRODUCT DESCRIPTION

The Mayo Abdominal Retractor is a foundational deep-cavity surgical instrument engineered specifically to provide broad, atraumatic soft-tissue and visceral displacement during major open abdominal and pelvic procedures. Designed with an expansive, smoothly curved paddle blade terminating in a soft retaining lip, this retractor allows surgical teams to hold back heavy muscle layers, adipose tissue, and viscera cleanly away from the operative corridor.

Built for lasting mechanical strength during major laparotomies, this dedicated Mayo Abdominal Retractor configuration features an overall tool length of 28 cm (\$109) paired with an expansive 70 mm blade width. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to ensure absolute durability, corrosion resistance, and a lifetime protection guarantee against rust, its broad blade architecture distributes manual counter-traction evenly to prevent parenchymal organ trauma. The ergonomic handle delivers superior pulling leverage and hand stabilization throughout multi-hour procedures.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Mayo Abdominal Retractor to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing general, gynecological, and urological operating theaters worldwide with unyielding clinical durability and dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)