



RETRACTOR MAKER

PRODUCT BROCHURE OF

Weitlaner Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Weitlaner Retractor

SKU RM-WR-SERIES-1

\$94 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
4" to 4 ½" (10–11.5 cm) — 2×2 or 2×3 prongs	RM-WR-SERIES-1-1	\$94
5" to 5 ½" (13–14 cm) — 3×4 prongs	RM-WR-SERIES-1-2	\$99
6 ½" (16.5 cm) — 3×4 prongs	RM-WR-SERIES-1-3	\$104
8" (20 cm) — 3×4 prongs	RM-WR-SERIES-1-4	\$109
9 ½" to 10" (24–25.5 cm) — 5×6 prongs	RM-WR-SERIES-1-5	\$114

PRODUCT DESCRIPTION

The Weitlaner Retractor is an internationally acclaimed self-retaining surgical spreader engineered specifically to provide reliable, hands-free soft-tissue and wound margin separation across general, orthopedic, plastic, and neurological procedures. Designed with outward-curving multi-prong jaws governed by a positive mechanical ratchet lock, this instrument holds incision borders open under continuous tension, eliminating assistant hand fatigue while preserving an unobstructed operative corridor.

Built to clinical framework specifications, this Weitlaner Retractor line features an extensive sizing matrix across five primary longitudinal reaches: 4" to 4 ½" (10–11.5 cm, \$94) with 2×2 or 2×3 prongs for small incisions, 5" to 5 ½" (13–14 cm, \$99) with 3×4 prongs for medium wounds, 6 ½" (16.5 cm, \$104) with 3×4 prongs for general and spine surgery, 8" (20 cm, \$109) with 3×4 prongs for deep tissue retraction, and 9 ½" to 10" (24–25.5 cm, \$114) with 5×6 prongs for maximum depth and heavy muscular exposure. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its prongs are available in sharp or blunt profiles to match tissue sensitivities.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Weitlaner Retractor to withstand repetitive hospital autoclave cycles, providing operating rooms, trauma suites, and outpatient surgery centers worldwide with exceptional mechanical durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)