



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cloward-Style Blade Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cloward-Style Blade Retractor

SKU RM-CSBR-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
13 × 75 mm	RM-CSBR-1-1	\$84
18 × 75 mm	RM-CSBR-1-2	\$89
20 × 75 mm	RM-CSBR-1-3	\$94
23 × 75 mm	RM-CSBR-1-4	\$99

PRODUCT DESCRIPTION

The Cloward-Style Blade Retractor is an indispensable handheld spinal surgical retractor engineered specifically to provide secure, non-traumatic exposure of the anterior cervical spine during disc and fusion interventions. Designed to gently retract the longus colli muscles and displace the esophagus and carotid sheath laterally, this instrument establishes a safe visual corridor to the intervertebral disc spaces.

Built for variable patient depths, this Cloward-Style Blade Retractor configuration features a uniform 19 cm overall chassis length paired with four selectable blade widths across a fixed 75 mm depth: 13 mm (\$84), 18 mm (\$89), 20 mm (\$94), and 23 mm (\$99). Surgeons can choose between a With Lip configuration designed to prevent muscle creep over the blade margin, or a Without Lip design for smooth, low-profile insertion in confined tissue planes. Its ergonomic flat handle provides steady, fatigue-free manual leverage.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Cloward-Style Blade Retractor to strict clinical tolerances, ensuring high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for hospital operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)